

Increased blood flow
generates a
natural healing response



Wound Therapy

Providing increased blood circulation
to promote wound healing naturally
from the inside

Supported by several peer reviewed papers and the recipient of multiple awards including best new product or service in Long Term Care by the Ontario Long Term Care Association.



For Canadian Use Only

Innovative Technology

A paradigm shift in the management of wounds

The cause:

- Wounds, including venous, arterial and diabetic ulcers often have impaired blood flow^{1,2}
- Impaired calf muscle pump function increases venous stasis and venous hypertension, and can negatively impact the severity of venous ulcerations^{3,4,5,6}


The treatment:

- Improved blood circulation results in enhanced wound closure^{2,7}, a natural healing response
- The geko™ device increases venous, arterial and microcirculatory blood flow in the lower limb in patients with chronic venous insufficiency^{8,9} and intermittent claudication¹⁰. It also reduces edema^{11,12}, augments the calf muscle pump¹³ and maintains TCpO₂ – promoting conditions suitable for wound healing^{14,17}

What is geko™ wound therapy?

Self-contained and wearable, the geko™ device:

- Low electrical current is required for nerve stimulation (versus the much higher power required to stimulate muscles) resulting in a pain-free experience
- Stimulates the common peroneal nerve, it activates the extensor muscles and stretches the antagonistic flexor muscles, acting as a calf muscle pump¹³
- Increases superficial femoral venous volume flow by 100%, femoral arterial volume flow by 75%¹⁶ and microcirculatory flux to the skin on the dorsum of the foot¹⁷ and thigh¹⁸ by 400%
- Increases blood flow equal to 60% of that achieved when walking¹⁵, so wearing it for 6 hours per day may give a physiological benefit comparable to over 3 hours of walking
- Benefits patients with chronic venous insufficiency over time^{8,9}
- Is simple, easy to use, small and lightweight (just 10g) battery operated (no leads or wires), enabling the patient to be as mobile and independent as possible
- Is worn for 6 hours per day, 6 days per week
- Through gentle muscle contractions, provides feedback so that patients feel engaged in their care, and may have better adherence to treatment protocols




Speckle spectroscopy¹⁹ – evaluation of a venous leg ulcer

As an example, when activated, the geko™ device caused a 225% increase in flux ($p < 0.001$) in the wound bed and a 67% increase in flux ($p < 0.001$) surrounding the peri-wound area²⁰. Increases in flux corresponds to an increase in microcirculatory blood flow, which is clearly seen in the comparison below. This increase in blood flow results in an increase in red blood cells carrying oxygen and nutrients necessary for energy metabolism and healing.


Further evidence can be reviewed at: www.gekowound.ca



Baseline speckle flow pattern



After activation of the geko™ device



Benefits of the geko™ device

The geko™ device increases venous and arterial blood flow while reducing pain²⁰ in individuals with lower leg ulcers.

In addition, consider the geko™ device²¹:

- In the management of lower leg edema that is contributing to reported pain
- In the management of stalled, chronic lower leg wounds that are not progressing along the expected healing trajectory (or wounds that can be predicted to be slow in healing from the onset)
- In conjunction with compression or when compression cannot be tolerated
- For patients with fixed ankle joints, those who are bedridden or those who have limited mobility

Preliminary evidence in an Ontario Home Care setting evaluation suggests this may be a first line treatment in conjunction with traditional therapy²².


A recent consensus document from Wounds Canada articulated the benefits of low frequency nerve stimulation (geko™) in the management of venous leg ulcers²¹

The authors hypothesized that the management of such patients is complex. Many clinical and patient factors must be considered or influenced for successful patient management and healing. The pathway to success is shown in the adapted illustration below, demonstrating the areas where the geko™ device has an influence. The table **on the next page** provides further supporting evidence.

How does the geko™ device overcome venous leg ulceration?

The model below illustrates complex factors that contribute to chronic venous insufficiency (CVI) and venous leg ulceration. Various stages of the wound physiology such as edema, pain, venous hypertension and decreased arterial flow can be treated with the geko™ device, contributing to the healing of venous and other lower leg ulcers. **The geko™ device ends this cycle.**

The geko™ device has an effect on each of the areas circled in red



Original created by Dee O’Sullivan-Drombolis, modified with permission for illustrative purposes²¹.

Research Evidence

The geko™ device has been the subject of scientific rigor to demonstrate its ability to increase blood circulation. The body of evidence continues to grow, targeting clinical issues such as CVI, in the management of lower leg wounds.

Clinical Issue with CVI	Device Effect
Abnormal Calf Muscle Pump	
<ul style="list-style-type: none"> Non-healing venous ulcers correlate with impairment of the calf muscle pumps 55% of patients with CVI have Calf Muscle Pump Dysfunction related to altered gait, causing venous hypertension⁶ 	<ul style="list-style-type: none"> The device creates concentric contraction of the extensor muscles that cause dorsiflexion of the ankle joint and passive stretch of the calf flexor muscles. The passive motion of the flexor muscle acts as a calf muscle pump, which may enhance venous return by increasing intramuscular pressure¹³ This may be effective in reducing venous stasis and edema, influencing muscle oxygenation¹³ The results may indicate that the geko™ device effectively counteracts increases in muscle blood volume and deoxygenated hemoglobin during venous stasis¹³
Edema	
<ul style="list-style-type: none"> Dependent edema begins in the peri-malleolar region and ascends the leg in early stages of CVI, changes over time to become fibrotic and indurated with Lipodermatosclerosis (LDS) due to changes in the fibrinolytic system²³ It may develop into veno-lymphedema 	<ul style="list-style-type: none"> Case studies have shown that some patients with chronic and complex edema have had edema reductions with the geko™ device¹¹ In a trial of the geko™ device with individuals with CVI, leg swelling reduced by 16% (p<0.05) in patients with venous disease^{8,9} Many patients with chronic venous ulcers who were not able to tolerate ANY compression therapy, or only tolerate minimal 10-15 mm Hg compression, have been able to either start or increase their level of compression therapy, leading to further edema reduction²⁴ There is also a fibrinolytic effect²⁵ with the geko™ device which may reduce the fibrotic changes of LDS
Incompetent Venous Valves	
<ul style="list-style-type: none"> 84% of people with VLUs have superficial vein valve failure²⁶ Failure of the deep vein valves speeds venous disease²⁶ and increases the risk of venous ulcers Both cause venous reflux and venous hypertension²⁶ 	<ul style="list-style-type: none"> The geko™ device reduces venous refilling and venous volume seen in venous stasis due to the activation of the muscle pumps²⁷ It decreases the amount of “sludge” blood (erythrocytes seen as light gray in an ultrasound image) that is not effectively ejected forward through the valves with cardiac systole/ diastole²⁸ When off, the Venous Sludge Index (VSI) was 53.5, when activated, the geko™ stimulation reduced the VSI to 7.6 (p=0.0005)²⁸
Decreased Range of Motion, Decreased Muscle Strength and Activation	
<ul style="list-style-type: none"> Decreased range of motion of the ankle can be related to nociceptive and neuropathic pain, woody fibrosis/LDS, edema, and fixed ankle joint related to CVI, over time develop decreased muscle strength and activation, and therefore decreased mobility^{5,29} 	<ul style="list-style-type: none"> In case series studies, patients have reported an increased ability to flex and dorsiflex their foot and ankle³³, with increased strength in their legs with increased exercise tolerance¹¹
Pain	
<ul style="list-style-type: none"> People living with VLUs often report pain as 10/10 and are unable to tolerate compression therapy, which is one of the key interventions in treating CVI²⁴ 	<ul style="list-style-type: none"> Up to 90% of individuals with chronic long-standing VLUs using the geko™ device indicated a marked reduction in pain and a subsequent reduction in narcotic usage^{20,30} Patients who were unable to tolerate therapeutic levels of compression pre-geko™ were able to start and/or increase to therapeutic levels with the effect of the geko™ device²⁴
Neuropathy	
<ul style="list-style-type: none"> Neuropathy in individuals with CVI without Diabetes is related to perineural degeneration, edema, and collagen replacement and contributes to trophic skin changes and impaired healing³¹ Worse in proximal medial and lateral malleolus, proximal medial and lateral calf and thigh³¹ Maybe an unrecognized source of pain in this population³¹ 	<ul style="list-style-type: none"> A study of a low frequency stimulation device to either the common peroneal or saphenous nerve³², depending on proximity to the ulcer, in conjunction with a four-layer compression bandaging system over 12 weeks, showed nearly four times greater improvement in the nerve sensation and two times the response to capsaicin applied topically, (both parameters reflecting improvement in C-fiber function)³² The improvement of C-fibre activation is also an indicator of the reversal of the neuropathy²¹
Decreased Arterial Flow	
<ul style="list-style-type: none"> 15 to 30% of people with venous disease will also have peripheral arterial disease (PAD)³³ 	<ul style="list-style-type: none"> The geko™ device augments arterial, venous and microcirculatory volume flow in PAD patients and may prove a useful treatment adjunct in these patients^{10,25}
Ambulatory Venous Hypertension	
<ul style="list-style-type: none"> Unabated venous hypertension may result in dermal changes with hyperpigmentation; subcutaneous tissue fibrosis, termed “lipodermatosclerosis”; and eventual ulceration²⁶ 	<ul style="list-style-type: none"> The geko™ device was tested in 19 healthy volunteers using settings of 100 μs, 200 μs and 400 μs while volunteers were standing, sitting and lying. Mean Venous Transit Times (VTT) from the dorsal foot to the popliteal vein were measured along with ambulatory venous pressure and leg volume The geko™ had a statistically significant impact, reducing VTT by up to 64%, Mean ambulatory pressure by up to 67% and leg volume by 17% (P< 0.001)²⁷

Clinical evidence – evaluation of the geko™ device in the management of venous leg ulcers

Venous leg ulcer^a

A 77 year-old female had a 10-year history of chronic ulcerations on the medial aspects of the lower legs, each episode lasting several months. She was concordant with high compression and had dressing changes 3 times a week. The geko™ device was used 5 days/week, 6 hours/day with time off. Pain reduced from 6/10 to 0/10 and the chronic wound gradually closed and she was fitted with compression stockings.



Prior to treatment – 10 year history



Closed at 10 months, 5 days

Painful leg ulcer^b

A 41 year-old female with a BMI >33kg/m², spontaneous leg ulcers, 6 weeks prior; required IV and later oral antibiotics; still on oral x 5 days at baseline. ABPI: Left 1.0, Right 1.2; Pain 10/10 initially. As wounds closed she graduated from low to high compression as pain decreased to 0/10. She was fitted with compression stockings.



Prior to treatment – 6 week history



Closed at 18 weeks

Non-healing venous leg ulcer^a

An 80 year-old female reported a 6.5 month history of VLU to the right medial malleolus, left medial malleolus and a pressure ulcer on her left heel. She was unable to tolerate compression due to pain and received nursing visits for wound care 3 times a week. geko™ therapy was increased from 2 hours/day to 4 hours/day. One wound closed in 18 days and the remaining in 2 ½ months. When pain was reduced she was fitted with compression. Her nurse commented on a change in her overall appearance and well-being.



Prior to treatment – 6.5 month history



Closed at 18 days

Unable to tolerate compression^a

A 74 year-old male with 10-year history of left leg VLU closures and recurrences with several infections, bilateral knee replacements and a history of a severe burn to the area as a child. Pain was reported as 8/10 at baseline, he could not tolerate compression. But with geko™ pain reduced to 5/10 and he graduated from low to more moderate compression. He paused geko™ treatment periodically to manage dermatitis. The wound closed following 9 months of treatment and he was fitted with compression stockings.



Prior to treatment – 10 year history



Closed at 9 months

Clinical evidence – evaluation of the geko™ device in the management of other types of wounds

Non-healing surgical amputation^b

A 77 year-old male, with CVI and diabetes, had a non-healing surgical amputation site of one toe on the right foot, 4.5 months in age, with previous bypass surgery to this leg 7 years prior. Angioplasty was performed 1 month before amputation of the toe. He also had a venous ulcer on the right shin, which had doubled in size over 3 months. He was in an inelastic Unna's paste boot dressing. His nursing visits went from every 2 days at baseline, to every 3 days by week 3. Both wounds closed at 5 weeks.



Prior to treatment –
4.5 month history



Closed at 5 weeks

Diabetic foot wound^c

A female with type 2 diabetes and a non-healing second toe amputation; wedge resection and multiple plantar DFU which did not heal, following 1 year of wound care. She had 3+ peripheral edema below the knee. The geko™ therapy was 6 hours/day, 5 days/week. Edema reduced after 2 weeks and all wounds were closed following 4 weeks of geko™ treatment.



Prior to treatment –
1 year history



Closed at 4 weeks

Woody Fibrosis^b

A 67 year-old male with type 1 diabetes a several year history of bilateral VLU and recurrent blisters. He had three hospitalizations in the year prior to the geko™ with cellulitis and sepsis from his legs requiring IV Antibiotics. geko™ was initiated for 6 hours/day 5 days/week and within 2 weeks his legs were getting softer and he had increased ankle mobility. The recurrent blisters decreased in frequency and duration. During the time of the evaluation he experienced only 1 course of oral antibiotics and no hospitalization.



Patient using the geko™ device
– 14 month history



Some areas closed at 12 weeks;
rest closed by 9 months

Pressure injury^d

A 92 year-old female with Atrial fib, type 2 diabetes, benign hypertension, arthritis, glaucoma, and dementia. Wound etiology appeared to be pressure-related, and offloading and repositioning schedule in place. ABPI not available; suspected some arterial compromise. Right heel 0.9 x 0.6 cm, covered with scab. Left heel 2.1 x 1.7 cm, covered in eschar and dry scab, surrounded by hyperkeratotic skin. Wound duration of 4 months, healed in 3 months with the geko™ in combination with conservative sharp wound debridement and best practices.



Prior to treatment –
4 month history



Closed at 3 months

References

1. Clarke-Moloney M, Lyons GM, Burke PE, O'Keeffe D, Grace PA. A review of technological approaches to venous ulceration. *Crit Rev Biomed Eng.* 2005;33(6):511-56
2. NHS Choices Venous Leg Ulcers March 2015. <http://www.nhs.uk/Conditions/Leg-ulcer-venous/Pages/Introduction.aspx>
3. Milic DJ, Zivic SS, Bogdanovic DC, Karanovic ND, Golubovic ZV. Risk factors related to the failure of venous leg ulcers to heal with compression treatment. *JVasc Surg.* 2009;49(5):1242-7. Epub 2009 Feb 23. [http://www.jvascsurg.org/article/S0741-5214\(08\)02007-7/fulltext](http://www.jvascsurg.org/article/S0741-5214(08)02007-7/fulltext)
4. Araki CT, Back TL, Padberg FT, Thompson PN, Jamil Z, Lee BC, Duran WN, Hobson RW 2nd. The significance of calf muscle pump function in venous ulceration. *JVasc Surg.* 1994 Dec; 20 (6):872-7; discussion 878-9. [http://www.jvascsurg.org/article/0741-5214\(94\)90223-2/pdf](http://www.jvascsurg.org/article/0741-5214(94)90223-2/pdf)
5. O'Brien JA, Edwards HE, Finlayson KJ & Kerr G. Understanding the relationships between the calf muscle pump, ankle range of motion and healing for adults with venous leg ulcers: A review of the literature [online]. *Wound Practice & Research: Journal of the Australian Wound Management Association*, Vol. 20, No. 2, Jun 2012: 80-85
6. Williams KJ, Ayekoloye O, Moore HM, Davies AH. The calf muscle pump revisited. *J Vasc Surg Venous Lymphat Disord.* 2014;2(3):329-34. doi: 10.1016/j.jvsv.2013.10.053. Epub 2014 Jan 28
7. Mosti G, Iabichella ML, Partsch H. Compression therapy in mixed ulcers increases venous output and arterial perfusion. *JVasc Surg.* 2012 Jan; 55(1):122-8. [http://www.jvascsurg.org/article/S0741-5214\(11\)01816-7/fulltext](http://www.jvascsurg.org/article/S0741-5214(11)01816-7/fulltext)
8. Williams KJ, Babber A., Ravikummar R, Ellis M, Davies AH. Pilot Trial of neuromuscular stimulation in the management of chronic venous disease. 2 Posters from VEINS Conference, UK. 2014
9. Williams KJ, Davies AH. Pilot trial of neuromuscular stimulation in the management of chronic venous disease. *British Journal of Surgery.* 2015;102:20
10. Barnes R, Shahin Y, Tucker AT, Chetter IC. Haemodynamic efficacy of the geko™ electrical neuromuscular stimulation device in claudicants. Oral presentation at Society of Academic & Research Surgery, 2014 Annual Meeting (January 8/9, 2014), Cambridge University, England. http://www.surgicalresearch.org.uk/wp-content/uploads/2013/10/1A_Vascular_Surgery_1.pdf
11. Ingves MV, Power AH. Two Cases of Transcutaneous Electrical Nerve Stimulation of the Common Peroneal Nerve Successfully Treating Refractory, Multifactorial Leg Edema. *Journal of Investigative Medicine High Impact Case Reports.* October-December 2014; 1-4. Available at: <http://journals.sagepub.com/doi/abs/10.1177/2324709614559839>
12. Wainwright TW, Immins T, Middleton RG. Poster Physiotherapy UK, October 2014, Birmingham
13. Zhang Q, Styf J, Ekström L, Holm AK. Effects of electrical nerve stimulation on force generation, oxygenation and blood volume in muscles of the immobilized human leg. *Scand J Clin Lab Invest.* 2014 Aug;74(5):369-77
14. Warwick D, et al. Microcirculation in the foot is augmented by neuromuscular stimulation via the common peroneal nerve in different lower limb postures: a potential treatment for leg ulcers. *International Angiology* 2015 April; 34(2):158-65
15. Tucker AT, Maass A, Bain DS, et al. Augmentation of venous, arterial and microvascular blood supply in the leg by isometric neuromuscular stimulation via the peroneal nerve. *Int J Angiol* 2010;19:e31-e37
16. Williams KJ, Moore HM, M Ellis and Davies AH. Haemodynamic changes with the use of a neuromuscular stimulation device compared to intermittent pneumatic compression. *Phlebology.* Online 10 April 2014. <http://phl.sagepub.com/content/early/2014/04/10/0268355514531255>
17. Jawad H, Bain DS, Dawson H, Crawford K, Johnston A, Tucker AT. The effectiveness of a novel neuromuscular electrostimulation method versus intermittent pneumatic compression in enhancing lower limb blood flow. *J Vasc Surg: Venous Lymphat Disord.* 2014;2(2):160-5
18. Bahadori S, Immins T, Wainwright TW. The effect of calf neuromuscular electrical stimulation and intermittent pneumatic compression on thigh microcirculation. *Microvascular Research* 2017; 111: 37-41
19. Qin J, Shi L, Wang H, Reif R, Wang RK. Functional evaluation of hemodynamic response during neural activation using optical microangiography integrated with dual-wavelength laser speckle imaging. *J Biomed Opt.* 2014;19(2):026013-1 to-4
20. Harding KG. A New Innovation in Wound Treatment. Presenting at CAWD Conference 2016.
21. Orsted HL, O'Sullivan-Drombolis D, Haley J, LeBlanc K, Parsons L. The effects of low frequency nerve stimulation to support the healing of venous leg ulcers. *Canadian Association of Wound Care Consensus Paper – November 2016*
22. Waterloo Wellington Community Care Access Centre Pilot Evaluation. 2017. Data on file.
23. Vivas A, Lev-Tov H, Kirsner RS. In the Clinic: Venous Leg Ulcers. *Ann Intern Med.* 2016;165(3):ITC17-ITC32
24. Harris C, Duong R, Vanderheyden G, Byrnes B, Cattryse R, Orr A, Keast, D. Evaluation of a Neuromuscular Electrical Stimulation Device for Non-Healing Venous Leg Ulcers. *Int. Wound J.* 2017, In press
25. Barnes R, Madden LA, Chetter IC. Fibrinolytic effects of peroneal nerve stimulation in patients with lower limb vascular disease. *Blood Coagulation and Fibrinolysis.* 2016; 27:275-280
26. Eberhardt RT, Raffetto JD. Chronic venous insufficiency. *Contemporary reviews in cardiovascular medicine circulation.* 2014;130:333-346. <http://circ.ahajournals.org/content/130/4/333>
27. Khanbhai M, Hansrani V, Sultan J, Burke J, McCollum CN. The effect of neuromuscular electrostimulation on lower limb venous physiology. Academic Surgery Unit, Institute of Cardiovascular Sciences, Manchester Academic Health Science Centre. Society of Academic & Research Surgery, VASCULAR 1 Wednesday 7 January, 2015 09.45-11.20 <http://www.surgicalresearch.org.uk/sars-2015/>
28. Lattimer C, Azzam M, Kalodiki E. Common peroneal nerve stimulation reduces blood sludging in the popliteal vein standing and lying. 2016. http://www.gekodevices.com/media/128135/acp_2016_geko.pdf
29. Back TL, Padberg FT, Jr., Araki CT, Thompson PN & Hobson RW, 2nd. Limited range of motion is a significant factor in venous ulceration. *J Vasc Surg* 1995; 22(5):519-23
30. Brooke J, Loney A. The geko™ a neuromuscular electrostimulation (NEMS) device and its healing effect on diabetic foot and venous leg ulcers. CAET/WOCN Conference, May 2015
31. Padberg FT, Maniker AH, Carmel G, Pappas PJ, Silva MB, and Hobson RW. Sensory impairment: A feature of chronic venous insufficiency. *J Vasc Surg* 1999;30:836-43
32. Ogrin R, Darzins P, Khalil Z. The use of sensory nerve stimulation and compression bandaging to improve sensory nerve function and healing of chronic venous leg ulcers. *Current Aging Science.* 2009;2(1):72-80
33. Marston W. Mixed Arterial and Venous Ulcers. *Wounds.* 2011;23(12):351-356

Case Study References

- a. Harris, C et al. Refractory Venous Leg Ulcers: Observational Evaluation of Innovative New Technology, *Int. Wound J.* 2017, In press
- b. Harris, C et al. Evaluation of a Neuromuscular Electrical Stimulation Device for Non-Healing Venous Leg. *Int. Wound J.* 2017, In press
- c. Case Study from Perfuse Medtec Inc. Archives used with patient permission
- d. Harris C. et al. Using a Muscle Pump Activator Device for Non-healing Lower Leg Wounds in Long Term Care Residents. Manuscript for submission to *Int. Wound J.* 2017, In press

Distributed in Canada by:



Information / Demonstrations / Ordering:

Email: orders@perfusemedtec.com

For more information on the geko™ device and the clinical evidence, please visit our website: www.gekowound.ca

MWBRCAN0350V1

© Firstkind and Perfuse 2017

geko™ is a registered trade mark of Sky Medical Technology Limited
OnPulse™ is a registered trade mark of Sky Medical Technology Limited